

Prompt gamma-ray imaging for particle beam therapy using nuclear emulsion

Toshiyuki Toshito, Mitsuhiro Kimura, Hiroyuki Ogino

Nagoya Proton Therapy Center

Yuta Shibamoto

Nagoya City University

Osamu Sato, Mitsuhiro Nakamura

Nagoya University

This work is supported by JSPS KAKENHI Grant Number 16K15349.


This work was carried out by the joint usage / research program of the
Institute of Materials and Systems for Sustainability (IMaSS), Nagoya University.



Nagoya Proton Therapy Center

Particle therapy

Depth-dose



Particle therapy is much more sensitive to the range uncertainty !

In Japan

Proton	13	21000	Pediatric
Carbon	5	15000	Bone and soft tissue

100 patients/year
radiotherapy

In the world proton 60, carbon 10



Nagoya Proton Therapy Center





Range uncertainty

- Calculation
 - CT-value to material conv. ~3%
- Machine
 - Beam energy ~1 mm
- Patient setup
- Internal motion
 - respiration
 - **Anatomical change**



Prompt gamma-ray imaging - principle

Gamma camera
SPECT(single photon
emission computed
tomography)




Moteabbed Phys. Med. Biol. 56 1063 (2011)




Nagoya Proton Therapy Center

Prompt gamma-ray imaging



– experimental history




Min Appl Phys Lett 89 183517 2006



Kurosawa Current Appl Phys 12 364 2012



Richter Radio Oncol 118 232 2016




Yamaguchi NIM A 833 199 2016



Prompt gamma-ray imaging using nuclear emulsion


- Detect $\gamma \rightarrow e^+e^-$ conversion below 10 MeV
 - Low energy extension of GRAINE
- Superior background rejection
- No pile-up, No dead-time
 - All events can be recorded
- Limited angular resolution




Experiments for proof-of-principle


- Dec. 2016, Jan. 2017 6.25 cm x 5 cm x 10 films
- Aug. 2017 12.5 cm x 10 cm x 10 films

Water Phantom


 2×10^{11}
200 MeV
Proton beam





Expected yields
 3.4×10^4 e+e-pair
 2.9×10^4 either e+ or e- is greater than 1 MeV
 1.5×10^4 both e+ and e- are greater than 1 MeV
 4×10^2 one is greater than 10 MeV and the other is greater than 1 MeV



Pilot analysis

4 cm x 4 cm in both sides of all emulsion films were scanned by using HTS.
 7×10^7 tracks / film were detected as micro tracks.
 5×10^6 tracks / film were reconstructed as base tracks.

Base tracks in the upmost film



Summary

- Prompt gamma-ray imaging is expected to be useful to visualize dose in the patient.
- Challenging is tracking of low energy electrons/positrons.
- We are testing prototype nuclear emulsion gamma-ray camera systems for proof-of-principle.